


How to Download and Install DLM 317

Overview

DLM 317 brings QoF v18. This will be delivered automatically via a DLM download as usual. However, practices using a local area network (LAN) can install DLM 317 manually if required:

1. First, you must ensure that all users are logged off Vision. From the Vision front screen go to **Login – Show Users** . There should be no users showing in the window.



2. Go to: http://www.inps4.co.uk/my_vision/vua/dlm/index.html
3. From the DLM 317 section, click on the **DLM 317 Download** folder and download the file to your desktop.
4. Find the file on your desktop, right click and select **Extract All**.




5. Choose a location to extract the files to.

- Once extracted, go to the chosen location and double click on the **Setup.exe** file.
- If prompted with an Open File Security Warning screen, press **Run**.



Open File Security Warning - Run

- The Vision Upgrade Installation wizard is displayed. Click **Next** to continue.




Vision Upgrade Installation - Next

- The files will start to install on your system. You are then prompted with "Did the Software Registration form print correctly?". Select **Yes**.

Note - If you have any problems, please call the helpline between the hours of **07:00 – 19:00** on **0845 351 1820** or **020 7501 7105**.

10. You are then prompted with "Vision 3.01 has been successfully installed. Click Finish to exit." Click **Finish**.



Finish to exit

11. The installation is now complete. Please see DLM 317 user guide or Clinical Audit on-screen help for changes.